Extremely precise. Diaper or Specimen Scale.

- Compact | Internal rechargeable battery Easy-to-clean | Three models
- Doran. Powering Performance.

doranscales.com



dorani

CUASS 111 10 × 0.005 15; 160 × 0.1 02: 4500 × 2 5; 4.5 × 0.002 kg

PC-400 DIAPER OR SPECIMEN SCALE

Compact, portable, easy-to-clean, and constructed entirely of stainless steel. Precision weigh in as fine as a 0.1 gram with performance features including a Quad-Spring shockproof base Auto Shut-Off, and fast battery recharging.

PC-400 SERIES

Diaper or Specimen Scale

The PC-400 is part of a complete line of high-quality Doran medical scales, designed to be durable and highly accurate, with innovative safety and performances features that support patient comfort and easy operation.

Constructed entirely of 304 stainless steel for easy wipe down cleaning. Extremely precise, with weight-in as fine as 0.1 gram. Large LCD read-out display. Power is provided by either an internal rechargeable battery pack, charged with the included AC adapter. The auto shut-off feature maximizes battery life.

INNOVATIVE DESIGN AND PERFORMANCE FEATURES

- Weight is quickly locked-on and displayed with Doran exclusive Weight Lock technology.
- Quad-Spring shockproof base design that virtually eliminates damage from shock or overload.
- · Lightweight and portable with exceptional build quality.
- Programmable auto shut off feature.
- Battery or AC power for continuous operation. Low Battery indicator.
- · Adjustable feet to maintain stability.
- · Quickly switch between lb, kg, g, oz, or lb:oz.

SPECIFICATIONS

doran

| Battery: 90 Days |
|--|
| Internal rechargeable battery and wall mounted UL Approved AC transformer standard |
| 1000 Weighments |
| On/Off, Zero, and Units |
| lb, kg, g, oz, or lb:oz |
| 304 Stainless Steel |
| Overall : 8.75" L x 6.5" W x 2" Platter : 6.5" x 6.5" |
| Model PC400-02 2 x 0.0002 lb (1000 x 0.1 g) |
| Model PC400-05 5 x 0.002 lb (2300 x 1 g) |
| Model PC400-10 10 × 0.005 lb (4500 × 2 g) |
| Model PC400-02 +/- 0.0002 lb (0.1 g) |
| Model PC400-05 +/- 0.002 lb (1 g) |
| Model PC400-10 |
| |

THREE-BUTTON CONTROLS

POWER On/Off with adjustable automatic shutoff timer for battery power applications.

ZERO used to zero-out the weight of any element to be included in weight of the fluids, stool or specimen.

UNIT button to set weight by lb, kg, g, oz, or lb:oz.

