

Plastic Mechanical Height Rod

KEY FEATURES

Accuracy and ease-of-use are the hallmarks of Doran’s medical scale line – and the DS1100 Mechanical Height Rod is no exception. Simply slide the height measurement arm into position. The patient’s height is shown in a large window just below the measurement arm. The bright red line makes the measurement very accurate and easy to read.

Precisely marked in centimeters and inches, the DS1100 stadiometer can be used to measure adult or pediatric patients.

The DS1100 is also very easy to install. The height rod mounts on the wall using three screws which are concealed by the durable plastic measuring tape.



Doran’s height rod features a user-adjustable measuring range, allowing customers to set the height during installation. The measuring tape is mounted so the 0 point aligns with the floor, and tape clips secure the height rod at the top and bottom. This flexible setup ensures accurate, easy-to-read measurements tailored to your specific height range needs.

SPECIFICATIONS

Warranty:	Two Years
Measuring Range:	Minimum Range 1.25 – 55.0" (3.0 – 140.0 cm) Maximum Range 38.5 – 92.5" (98.5 – 235.0 cm) Typical Range 28.0 – 82.0" (71.0 – 208.0 cm)
Graduations:	1/16" (1 mm)
Construction:	ABS Plastic
Dimensions:	Height Rod 59" L x 3" W x 1" D Measurement Arm 11" L x 8" W x 3.5" D

